



Wound Closure Strips

Description

Can be used in conjunction with or as an excellent alternative to sutures or staples

Features & Benefits

- Adhesive backed, thin non-woven, reinforced fibers that are arranged in a bilateral pattern for high tensile strength
- Provides support following early suture/staple removal
- Keeps wound closed without significant tension, edema, stretch or movement
- Translucent and porous
- Quick, simple application
- Non-invasive, less physical trauma
- Comfortable to wear
- Improved esthetic results; no puncture marks

Indications:

- Provide support following early suture/staple removal
- Closed wounds without significant tension, edema, stretch or movement
- Secure devices to skin (such as IV)
- Secure skin grafts in burn patients
- Secure skin flaps

Contraindications:

- Infected wounds



Connect with us!

Find our full line of
products online at:
www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779





PASSION. PARTNERSHIP. POSSIBILITIES.

Item	Description	Packaging
Sterile - Latex Free		
5150	1/8" x 3"	5/pk, 50 pk/bx, 4 bx/cs
5151	1/4" x 3"	3/pk, 50 pk/bx, 4 bx/cs
5152	1/4" x 1 1/2 "	6/pk, 50 pk/bx, 4 bx/cs
5156	1/4" x 4"	10/pk, 50 pk/bx, 4 bx/cs
5157	1/2" x 4"	6/pk, 50 pk/bx, 4 bx/cs
5158	1" x 4"	4/pk, 50 pk/bx, 4 bx/cs
Item	Description	Packaging
Non-Sterile - Latex Free		
C-02	1/4" x 3"	500/bx
C-03	1/4" x 1 1/2 "	500/bx
C-04	1/4" x 4"	500/bx
C-05	1/2" x 4"	500/bx
C-06	1" x 4"	500/bx

Connect with us!

Find our full line of
products online at:

www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation

2 Fleetwood Court,

Ronkonkoma, NY 11779

