



ELASTIC BANDAGES

Elastic bandages come in a variety of lengths, sizes and material options; as well as closure options such as clips or hook & loop. The intended design of an elastic bandage allows the bandage to stretch which causes localized pressure. The strength and compression features will vary from product to product to meet the needs of the care provider.



Premium
(Hook & Loop)



Standard
(with Clips)

Elastic bandages, also known as compression bandages, are often used to put gentle pressure on the tissue around an injury, as well as providing additional support to the injured area.

Other uses for elastic bandages are:

- To help decrease pain and reduce swelling of an injured area
- To hold wound dressings in place
- To wrap around the splinted arm or leg during healing
- To improve blood flow to a limb
- To hold cold or hot compresses in place

Best practices when using elastic bandages:

- Overlap the elastic bandage by one-half to one-third of its width each time you go around
- Do not wrap the elastic bandage too tight as to negatively impact the flow of blood
- To address good blood flow you may want to take off the wrap periodically for a few minutes provided it is in accordance with your medical caregiver instructions.
- If you experience numbness, tingling associated with use of the bandage or the area feels cold and has a blue tint it may be an indication of it being wrapped too tightly. Unwrap area, rub gently and when conditions go back to normal re-wrap the area.

See other side for sizes & styles available.

Quality.

Value.

Performance.



Premium Elastic Bandages

- Designed with hook and loop closure at both ends; eliminates the need for clips
- The crotched construction provides a better hook to secure and minimizes pulls on the fabric for repeat use
- Ideal for high compression situations
- Individually wrapped
- Perforated carton for easy dispensing

Non-Sterile - Premium Elastic Bandages - Self Closure

ITEM #	DESCRIPTION	PACKAGING
512	2" x 5yds	1/pk, 10/bx, 5 bx/cs
513	3" x 5yds	1/pk, 10/bx, 5 bx/cs
514	4" x 5yds	1/pk, 10/bx, 5 bx/cs
516	6" x 5yds	1/pk, 10/bx, 5 bx/cs
512LF	2" x 5yds, Latex Free	1/pk, 10/bx, 5 bx/cs
513LF	3" x 5yds, Latex Free	1/pk, 10/bx, 5 bx/cs
514LF	4" x 5yds, Latex Free	1/pk, 10/bx, 5 bx/cs
516LF	6" x 5yds, Latex Free	1/pk, 10/bx, 5 bx/cs

Sterile - Premium Elastic Bandages - Self Closure

ITEM #	DESCRIPTION	PACKAGING
522LF	2" x 5yds, Latex Free	1/pk, 50 pk/cs
523LF	3" x 5yds, Latex Free	1/pk, 50 pk/cs
524LF	4" x 5yds, Latex Free	1/pk, 50 pk/cs
526LF	6" x 5yds, Latex Free	1/pk, 50 pk/cs

Standard Elastic Bandages

- Knitted construction offers strength, comfort, and excellent breathability
- Provides firm compression and consistent support
- Individually wrapped and includes two metal clips
- Perforated carton for easy dispensing

Non-Sterile -Elastic Bandages - with Clips

ITEM #	DESCRIPTION	PACKAGING
502	2" x 4.5yds	1/pk, 10 pk/bx, 5bx/cs
503	3" x 4.5yds	1/pk, 10 pk/bx, 5bx/cs
504	4" x 4.5yds	1/pk, 10 pk/bx, 5bx/cs
506	6" x 4.5yds	1/pk, 10 pk/bx, 5bx/cs
502LF	2" x 4.5yds, Latex Free	1/pk, 10 pk/bx, 5bx/cs
503LF	3" x 4.5yds, Latex Free	1/pk, 10 pk/bx, 5bx/cs
504LF	4" x 4.5yds, Latex Free	1/pk, 10 pk/bx, 5bx/cs
506LF	6" x 4.5yds, Latex Free	1/pk, 10 pk/bx, 5bx/cs