

Passport® 12

PATIENT MONITOR

Features and Benefits

- 12" high resolution touchscreen display
- Up to 8 waveforms
- Standard features include 3 or 5-lead ECG, Masimo SET® SpO₂, 2 IBP, NIBP, respiration, temperature, integrated recorder
- Optional Nellcor® OxiMax® SpO₂ available
- Additional module parameter options include:
 - Sidestream CO₂
 - 2 IBP
 - 2 IBP + CO
 - 2 IBP + CO + Sidestream CO₂ or Microstream® CO₂
 - Multi-gas + 2 IBP + CO
- Available ST and arrhythmia analysis
- Extensive data storage capability for trend data, alarms, events, NIBP measurements and up to 48 hours of full disclosure
- Available wireless 2.4 GHz (802.11b/g/n)
- Connects to DPM Central Station or Panorama® Patient Monitoring Network
- Integrated HL7 and DIAP protocol interfaces
- Compatible with Mindray eGateway for direct HL7 communication to EMR
- Defibrillator synchronization
- Ergonomic carrying handle and multiple mounting options offer unprecedented mobility



your powerful partner across the care continuum

Mindray proudly offers the Passport 12, an inspired result of ingenuity and efficiency offering a modular approach to the traditional configured multi-parameter monitor. Flexible for virtually any clinical environment, the Passport 12 offers a plug-and-play modular structure allowing for more advanced monitoring functions such as CO₂, two additional invasive pressures, cardiac output, or a combination of all three. Create a compact, yet powerful OR monitoring solution when utilizing the available multi-gas module that automatically identifies and measures the inspired and expired values of O₂, CO₂, N₂O and 5 anesthetic agents.

Passport 12 extends beyond just the device itself by offering integrated HL7 and DIAP communication protocols. Wired and wireless networking capabilities can also be integrated, enabling connectivity to the DPM Central Station and Panorama Patient Monitoring System for a centralized patient data system.

Delivering enhanced capability, more optional features and extraordinary value, the Passport 12 fits your acute-care monitoring needs.

mindray
NORTH AMERICA

Passport® 12

Flexible Modular Design and Comprehensive Monitoring

The Passport 12 provides a modular approach to the traditional configured multi-parameter monitor. With a plug-and-play modular structure, the Passport 12 offers advanced monitoring functions such as IBP and CO₂, adding flexibility according to your clinical requirements. Additionally, these modules can be shared with the Passport 8 monitor, offering an expandable clinical solution at an exceptional value.



2 IBP



Sidestream CO₂



Sidestream CO₂ + 2 IBP + CO



Multi-gas Module + 2 IBP + CO



Integrated bed rail hook



Plug-and-play modules



Passport 8 and 12 with View Other

Interoperability for Improved Workflow

- Centralize patient data using a wired or wireless connection to the DPM Central Station and the Panorama Patient Monitoring Network
- View Other feature allows real time patient data to be viewed between the Passport 8 and Passport 12 bedside monitors. Up to 10 beds can be viewed simultaneously, whether connected to a central station or not
- When networked with the DPM Central Station, the Passport 12 can view telemetry devices at bedside
- Integrated DIAP and HL7 communication protocols for data export to EMR

Convenient Ergonomic Design

- Compact, ergonomic design allows modules to be easily inserted and released – with just one hand
- High-resolution LCD screen with intuitive touchscreen
- Customizable quick keys provide quick, convenient access to frequently used functions
- Integrated ECG, Resp, IBP, SpO₂, NIBP, Temp and 3-channel recorder
- Plug-and-play modules increase flexibility in monitoring
- 360-degree visible alarm indicator
- Convenient bed rail hook for ease of use during transport

Mindray DS USA, Inc.
800 MacArthur Blvd., Mahwah, NJ 07430
Tel: 800.288.2121 Tel: 201.995.8000 Fax: 800.926.4275 www.mindray.com

©2015 Mindray DS USA, Inc.
Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. All brands and product names are trademarks of their respective owners. Specifications subject to change without notice. Contact your Mindray Representative for the most current information. 05/15 P/N: 0002-08-8006 Rev D

mindray
NORTH AMERICA