

SallyTube Single Patient Use

Go green and transfer patients and residents safely with SallyTube recyclable slide tubes.

SallyTube® *Single Patient Use*


SallyTube Single Patient Use slide tubes can greatly improve infection control as well as improve manual handling over reusable slide sheets. With generous 50" x 48" standard dimensions, the SallyTube can accommodate patients and residents of all sizes.

During general patient or resident use, SallyTube Single Patient Use slide tubes can be used for several transfers and simply wiped down as needed or disposed of after use in high infection areas or isolation.




Standard Sizes

ST-5848S (58" x 48")



ST-5829S (58" x 29")



SallyTube Single Patient Use slide tubes are the perfect solution for lateral transfers in radiology, emergency departments, operating rooms, patient transport and isolation.



Ideal for lateral transfers of patients from to bed to bed, or stretcher to stretcher.



An optional plastic dispenser can be mounted in each room for convenient access.

SallyTube Single Patient Use



STEP 1

Pull out one sealed SallyTube from easy access wall dispenser



STEP 2

Insert SallyTube under patient

STEP 3

Transfer patient to stretcher / bed / examination table




Illustration above demonstrates use of SallyTube® (ST-5848S) with Tollos™ Draw Sheet with Handles (DPP7533H-80)

SallyTube Single Patient Use

- SallyTube single patient use lateral transfer slide tubes are the perfect solution for lateral transfers in radiology, emergency departments, operating rooms, patient transport and isolation.
- Caregivers can quickly slip the SallyTube under the patient or resident for repositioning, pivoting and lateral transfers.
- Made from a recyclable material, SallyTube single patient use slide tubes can be disposed of after use.

Plastic Dispenser for Convenient Access

An optional plastic dispenser can be mounted in each room for convenient access.

Features

- Improved handling for lateral transfers and repositioning
- Heat sealed edges for extra strength and infection control
- No laundering cost
- Generous 58" x 48" dimensions
- Made from recyclable material
- Also available in reusable material

TECHNICAL SPECIFICATIONS

Length (Overall)	58" (145cm)
Width (Overall):	48" (125cm)
Weight Capacity:	794lbs (360kg)
Recyclable:	Yes
Package:	Case of 5 boxes, Box of 10
Custom Sizes:	Upon request
Standard Sizes (by Model Number):	
ST-5848S	58" (145cm) x 48" (125cm)
ST-5829S	58" (145cm) x 29" (73cm)

www.tollos.com

U.S. Corporate Office
1 Easter Court, Suite J
Owings Mills, MD 21117
410.363.1515 PHONE
888.220.4566 TOLL FREE

Canadian Office
75 Dymont Road
Barrie Ontario L4N 3H6
705.733.0022 PHONE
800.565.7075 TOLL FREE