

arlington

recliner

rc1294


dimensions

29w . 36d . 44h



healthcare furnishings

|durfold|



arlington recliner

rc1294

front and rear swivel casters for ease of steering and mobility

3" locking rear casters for fixed positioning

seat height positioned to allow for ease of sitting and standing for elderly users and patients

ctb-133 tested for compliance

footpiece release is weight activated thereby preventing bending to mechanism from excessive and improperly placed weight on footpiece

wall-saver design to allow more open floor space for traffic and circulation

three position mechanism

prominent head pillow for support in a variety of resting positions

* durfold reserves the right to make changes to products or information at any time without notice

	recliner
w	29
d	36
h	44
sw	22
sd	23
sh	21
ah	28
rd	69

* rocker recliner available
nominal dimensions given

|durfold|

durfold corporation
102 upton drive, jackson, ms 39209
p.800.345.6849 f.601.922.6244
durfold.com