

# QUICKVUE<sup>®</sup>+

One-Step  
hCG Combo **TEST**

## A Simple Test You Can Depend On For The Early Detection Of Pregnancy



### Rapid Diagnosis

- Results in 5 minutes (serum), 3 minutes (urine).
- Positive results may become visible sooner.

### Excellent Performance

- >99% sensitivity; specificity; accuracy.
- 10 mIU/mL sensitivity in serum or 20 mIU/mL sensitivity in urine.

### Convenient

- Room temperature storage.
- 21-month shelf life.

### Easy to Use

- Simple, one-step procedure uses only 4 drops of sample.
- Two-color plus / minus (+ / -) endpoint for easy-to-read and interpret results.

### Reliable

- Highly sensitive assay with built-in controls.

### Lot-to-Lot Consistency


- Highly controlled and monitored manufacturing environment.

**QUIDEL<sup>®</sup>**

For more information on the QuickVue+ One-Step hCG Combo test kit, visit our website or call your Quidel Account Manager at 1.800.874.1517 (in the U.S.) or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [quidel.com](http://quidel.com)

## Procedure:



To order: Call your Quidel Account Manager or your local distributor.

QuickVue+ One-Step hCG Combo – 30 Test Kit:

Catalog # 00178

QuickVue+ One-Step hCG Combo – 90 Test Kit:

Catalog #00179

Serum hCG Control Set: Catalog #00281

QuickVue+ One-Step hCG Combo – 90 Test Kit:

Catalog #00179

hCG Control Set (Urine):

Catalog # 00272

# QUIDEL®

For more information on the QuickVue+ One-Step hCG Combo test kit, visit our website or call your Quidel Account Manager at 1.800.874.1517 (in the U.S.) or 858.552.1100 (outside U.S.).

**Quidel Corporation** • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [quidel.com](http://quidel.com)

Distributor Information: