

Technical Bulletin

QuickVue[®] Influenza and RSV Test Kits: Synagis[®] Testing

Synagis has been evaluated at a concentration of 1.0 mg/mL and was shown to not interfere with test performance of the following kits:

- QuickVue Influenza (00317)
- QuickVue Influenza A+B (20183)
- QuickVue RSV (20193)
- QuickVue RSV 10 (20222).

Please call Quidel Technical Support at 800-874-1517 (USA only), 858-552-1100 or technicalsupport@quidel.com if you have any questions regarding the QuickVue Influenza or RSV tests or any other Quidel product. Our hours of operation are Monday-Friday, 7:00 a.m.-5:00 p.m. Pacific Time.

You may also visit our website at www.quidel.com or www.flutest.com information on Quidel's line of Rapid Diagnostics, Bone Health and Autoimmune & Complement product lines. Other product information available on our website includes: CPT codes, CLSI, MSDS and package inserts.