

DISPOSABLE
RESTING ECG
ELECTRODE

TAB ELECTRODE

PROTAB

The BIOPROTECH PROTAB System
Your best choice for an Excellent Trace Quality in
Short Term Diagnostics

Today, facilities are challenged to reduce operating costs while maintaining high-quality patient care. The unique needs of resting ECG monitoring are met with the Bioprotech family of electrodes designed for this critical diagnostic application. Unique silver/silver chloride and conductive adhesive hydrogel of PROTAB provides stable baselines and a reliable signal tracing of consistent high quality.



PROTAB

Easy to apply and remove

The PROTAB provides quick and easy access to the electrode for single-handed application. The tab is convenient for leadwire placement and for ease of removal.

Time Efficiency

Unlike suction electrodes, PROTAB on a handy card are applied quickly on all types of patients. Strong performance gel requires no clean up after the ECG.

Superior adhesion

Excellent adhesion with the flexibility to securely reposition the electrode saves clinician time and facility cost.

Model	Backing type	Gel type	size	Radiolucent
PT-2334	Paper/PET	Conductive adhesive	23x34mm	Yes

**PT-2334**