

DISPOSABLE
ECG
ELECTRODE

TELECTRODE


How to Reduce Costs without compromising Quality?

Bio Protech has electrode solutions for every ECG application you need

Today, healthcare facilities are challenged to reduce operating costs while maintaining high-quality patient care diagnostics. Bio Protech can help you achieve this goal!

The Bio Protech TELECTRODE electrode

Bio Protech continually innovates new and improved solutions for ECG monitoring to meet the growing customer demands for affordable, easy to use, and accurate ECG electrodes. Bio Protech has developed a new range of uniquely designed ECG electrodes for all round diagnostic and monitoring applications whatever your specific requirements.




TELE CTRODE


Product Line

※ Model : T○○○C = Radiolucent , T○○○W = Wet Gel


Model	Size	Backing	Gel	Packaging
T715	55mm Round	PE Foam	Solid Hydrogel (Φ16)	50 pcs / Pouch 20 Pouchs / Box 10 Boxes / Case
T716	43X45mm Teardrop	PE Foam	Solid Hydrogel (Φ16)	
T815	55mm Round	Cloth	Solid Hydrogel (Φ16)	
T816	43X45mm Teardrop	Cloth	Solid Hydrogel (Φ16)	
T832C	32mm Round	Cloth	Sticky Full Gel	
T916	43X45mm Oval	Clear tape	Solid Hydrogel (Φ16)	




T715




T716




T815



T816



T815



T916C



T832C