

Reusable Purse String Clamp

Reusable, stainless steel, ring-handled instrument

Designed to assist in the mechanical construction of a purse string suture.

Reusable Purse String Clamp features

- Closes and secures tissue in a purse string pattern
- Can be used on any part of the gastrointestinal tract except the esophagus
- Five-year product life expectancy under normal usage conditions



Ordering Information

Product Code	Description	Qty/box
EH40	Reusable Purse String Clamp	1
D-9080	Sterile 2-0 prolene polypropylene suture, blue monofilament, 30 inch, 2 needles, straight taper point (malleable)	12

How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

Johnson & Johnson Gateway:
www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

Global Healthcare Exchange (GHX):
www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI):
JJHCS EDI Help Line: 1-800-262-2888

Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at **customersupport@eesus.jnj.com** or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.