

Curved to go deep. Efficient transection.

CONTOUR® Curved Cutter Stapler



Unique curved head

- Unique curved head designed to cut and staple deep in the pelvis
- Parallel closure and tissue retention pin designed to transect the rectum in a single perpendicular firing
- CONTOUR® allows placement of a 40mm cut line in the width of a 30mm space

More precise transection

- CONTOUR, with its simultaneous firing, may eliminate one reload per procedure¹
- 78% less tissue slippage during firing compared to Covidien Endo GIA™ Radial Reload with Tri-Staple™ Technology¹

Applicable specialties

- Colorectal
- General
- Thoracic

¹ Benchtop testing on porcine rectum. Normalized mean tissue movement along the cut line from before clamping on tissue to after firing the CONTOUR® Curved Cutter Stapler (CS40G) vs ENDO GIA™ ULTRA Handle (EGIAUSTND) and Endo GIA™ Radial Reload with Tri-Staple™ Technology (EGIARADMT) on 2.1mm to 2.9mm thick tissue. 3.8mm (9.4%) vs 17.4mm (32.5%), p<0.001.

Product overview

- Four curvilinear rows of staples with a single cut between
- Allows low pelvic access without handle obstruction over existing linear staplers
- Parallel closure for controlled compression and less tissue movement
- Tissue retaining pin with manual closure option holds tissue in place during compression and firing
- Breakthrough washer provides audible and tactile feedback every firing
- Devices are reloadable and may be fired up to 6 times in a single procedure



Device code	Tissue	Staple rows	Cut line length	Staple line length	Max firings	Qty/box
CS40B	■ Regular	4	40mm	51mm	6	3
CS40G	■ Thick	4	40mm	51mm	6	3



Reload code	Tissue	Closed staple height	Open staple height	Wire diameter	Qty/box
CR40B	■ Regular	1.5mm	3.5mm	0.20mm	6
CR40G	■ Thick	2.0mm	4.7mm	0.23mm	6

Ordering info

How to order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

Electronic ordering options

(Note: Placing orders electronically avoids minimum order fees for your hospital).

Johnson & Johnson Gateway: www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-886-JNJGATE (565-4283).

Global Healthcare Exchange (GHX): www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI): JJHCS EDI

Help Line: 1-800-262-2888

Nonelectronic/Manual ordering options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact Customer Support Center at customersupport@eesus.jnj.com or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.