

Thin. Nimble.  
And **incredibly**  
**precise.**



Introducing the  
**ECHELON FLEX™ Powered Vascular  
Stapler with Advanced  
Placement Tip**

Enables the **most precise placement**  
on fragile vessels

- **26% narrower curved, blunt tip anvil**, providing superior maneuverability around vascular structures
- **26% thinner shaft**<sup>1</sup> offering the greatest angle of reach in the tight intercostal space
- **11% greater manual articulation in each direction**,<sup>2</sup> allowing more flexibility during final placement

A narrower anvil, thinner shaft and more articulation than Endo GIA™ Curved Tip Reload with Tri-Staple™ Technology.

**Enhanced stability** on fragile vessels

- **83% reduction in tip movement**<sup>3</sup> during firing for less movement during transection, potentially resulting in less trauma

**Examples of applicable procedures**

Designed specifically for transection of fragile and/or difficult-to-reach vessels

- **Thoracic**—lobectomy, segmentectomy, pneumonectomy
- **Solid organ**—nephrectomy, donor nephrectomy, hepatectomy

1 Approach angles assessed in a virtual CAD environment in the 5th intercostal space.

2 Based on articulation data from IFUs of each device.

3 Benchtop testing on porcine carotid arteries. Surgeons (n=24) fired each instrument / reload once: PSE45A / ECR45W and EGAUSTND/ EGIA45AVM. Distal tip motion measurement during the firing cycle showed a median reduction of 83% and a range of reduction of 53% to 96% in tip movement of PSE45A/ ECR45W vs. EGAUSTND/ EGIA45AVM.

# Product components

- ECHELON FLEX™ Powered Vascular Stapler with Advanced Placement Tip
- ENDOPATH ECHELON™ Vascular White Reload for Advanced Placement Tip



## ECHELON FLEX™ Powered Vascular Stapler

CODE	DESCRIPTION	QUANTITY PER SALES UNIT
PVE35A	35mm Standard Articulating Endoscopic Linear Cutter,* 320mm	3

\*Instrument does not contain a reload.



## ENDOPATH ECHELON™ Vascular White Reload for Advanced Placement Tip (35mm, 4row)

CODE	COLOR	DESCRIPTION	OPEN STAPLE HEIGHT	STAPLE ROWS	CLOSED STAPLE HEIGHT	QUANTITY PER SALES UNIT
VASECR35	White	White reload intended for use in vascular/thin tissue	2.5mm	4	1.0mm	12

## How to Order

### Electronic Ordering Options

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS).

The following electronic order placement methods are preferred:

- J&J Gateway—[www.jnjgateway.com](http://www.jnjgateway.com) (for questions about your order, visit the Web site or call 1-866-JNJ-GATE)
- Electronic Data Interchange—JJHCS Help Line: 1-800-262-2888

### Nonelectronic/Manual Ordering Options

JJHCS—Call 1-800-255-2500 between the hours of 8:30 am and 8:00 pm Eastern time, or fax your order to 1-732-562-2212.

## Customer Support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter, or complaints, please contact our Customer Support Center by calling 877-ETHICON (384-4266). Our support center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.

Visit [www.ethicon.com/poweredvascular](http://www.ethicon.com/poweredvascular) for more information about the ECHELON FLEX Powered Vascular Stapler.

**For complete product details, see Instructions for Use available at [www.e-ifu.com](http://www.e-ifu.com).**