

Enseal®

ENSEAL® G2 **Curved and Straight** Tissue Sealers


Secure sealing and versatility
in the palm of your hand.



ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping
the future
of surgery

Secure **sealing** and **versatility** in the palm of your hand




Secure

- More than twice the compression at the distal tip of the jaw than LigaSure™ 5mm¹
- Less proximal end bleeding† and more consistent seal strength over multiple activations‡ than LigaSure™ 5mm
- Seals vessels up to 7mm diameter and lymphatics²

Versatile

- Delivers energy to tissue even with the jaws open, enabling spot coagulation, enhancing dissection and creating otomies

ENSEAL® G2 Tissue Sealers vs LigaSure™



† In a preclinical study utilizing an animate porcine mesentery model, ENSEAL® G2 Curved and Straight had significantly less bleeding at the proximal end than LigaSure™ 5mm (0.0% vs. 24.0%), P=0.001.

‡ In a benchtop study with ENSEAL® G2 Curved and Straight using 5mm porcine carotid arteries, vessel burst pressure was measured over 52 activation cycles P<0.05.

1. In side-by-side benchtop testing, average compression of distal 0.11 in of jaws used in ENSEAL® G2 Curved was 187 psi (standard deviation 81 psi) and ENSEAL® G2 Straight was 104 psi (standard deviation 36 psi) vs 42 psi (standard deviation 8 psi) for LigaSure™ 5mm. P<0.0000. 2. Data on file at Ethicon Endo-Surgery, Inc. Studies conducted on porcine mesentary.

Secure **sealing** and **versatility** to do more

A comprehensive redesign — in both the handle and activation sequence — ensures comfortable and confident operation.

Integrated activation


- Provides energy activation prior to cutting

Audible tones

- Provides feedback throughout the seal cycle

Ergonomic design

- Designed for comfort and performance



The Ethicon Endo-Surgery Generator
Revolutionary power
for the ENSEAL® G2 Tissue Sealers (sold separately)



ENSEAL® G2 Tissue Sealers

The complete ENSEAL® G2 Tissue Sealer Portfolio

All ENSEAL® and HARMONIC® devices are compatible with the Ethicon Endo-Surgery Generator.

ENSEAL® G2 Curved



Code	Shaft diameter	Shaft length	Qty
NSLG2C14	5mm	14cm	6
NSLG2C25	5mm	25cm	6
NSLG2C35	5mm	35cm	6
NSLG2C45	5mm	45cm	6

ENSEAL® G2 Straight



Code	Shaft diameter	Shaft length	Qty
NSLG2S14	5mm	14cm	6
NSLG2S25	5mm	25cm	6
NSLG2S35	5mm	35cm	6
NSLG2S45	5mm	45cm	6

ENSEAL® G2 Curved and Straight Tissue Sealers – Secure sealing and versatility to do more.

To experience the ENSEAL® G2 Tissue Sealers, contact your Ethicon Sales Professional
Call 1-877-ETHICON. Connect at www.ethicon.com/ensealg2

ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping
the future
of surgery

Enseal®