


Harmonic

Product Literature:
HARMONIC SYNERGY®
Combination Hook Blade



HARMONIC SYNERGY® Combination Hook Blade

The HARMONIC SYNERGY® Combination Hook Blade

HARMONIC® Technology uses ultrasonic vibration to simultaneously cut and coagulate soft tissue. This minimizes thermal tissue damage, may eliminate the need for electrocautery and may provide important recovery advantages to your patients.

Precise: Minimal lateral thermal tissue damage for safer dissection near vital structures

Reliable: Seals and divides vessels $\leq 2\text{mm}$, as well as lymphatics

Versatile: Cuts, coagulates and dissects, reducing instrument exchanges



HARMONIC SYNERGY® Combination Hook Blade

Product facts

Procedure Applications

Procedure applications include plastic surgery, exposure to orthopedic structures (such as spine and joint space), ENT and breast procedures, when bleeding control and minimal thermal injury are desired.

Plastic

Abdominoplasty, breast reduction

ENT

UPPP, tonsillectomy

Breast

Mastectomy, lumpectomy

Spine and Orthopedics


Used for soft tissue access to orthopedic structures



HARMONIC SYNERGY® Combination Hook Blade

Tissue effect comparison

Note: Tissue difference with HARMONIC® Technology



Bovie®

Tonsillectomy using electrosurgery



HARMONIC® Technology

Tonsillectomy using HARMONIC® Combination Hook Blade, HK105

HARMONIC SYNERGY® Combination Hook Blade

Specifications & product codes

The HARMONIC SYNERGY® Combination Hook Blade with Torque Wrench



Product Code	Quantity/Sales Unit	Price
SNGHK2	6	Contact Your Ethicon Sales Professional


HARMONIC SYNERGY® Combination Hook Blade

The HARMONIC SYNERGY® platform


One HARMONIC SYNERGY® Platform, Three Blade Options



Combination Hook Blade
SNGHK2



Hook Blade
SNGHK



Curved Blade
SNGCB