

Controlled compression from thin to thick tissue.

Proximate® Linear Cutter



Superior staple form in thick tissue compared to any other leading brand of open linear cutter¹

- 51% fewer malformed staples than Autosuture DST Series GIA (GIA8038S)

Reduced risk of leaks

- Proximate TLC (TLC55 and TLC75) with its unique CAM mechanism and tissue retention pin helps reduce staple line leakage compared to Autosuture DST Series GIA (GIA8038S)²

Applicable specialties

- Colorectal
- Gastric
- General
- Thoracic

¹ Benchtop testing in porcine tissue vs AutoSuture DST Series GIA (GIA8038S), Frankenman (LC80B), Reach (RLC8038S), Touchstone (LC8038) at 2mm-5mm thickness range (p<0.001).


² Benchtop testing of luminal leak pressure in an excised pre-clinical tissue model.

Product overview

- Unique CAM mechanism and tissue retention pin help ensure consistent staple formation
- Four different reloads available to accommodate a variety of tissue thicknesses
- Three lengths available to accommodate a wide range of procedural applications
- Intermediate locking position allows the surgeon to manipulate tissue without separating the device



Description	Device code	Tissue	Staple rows	Cut line length	Staple line length	Number of staples	Max firings	Qty/box
Proximate Linear Cutter 55mm	TVC55 TLC55 TCT55	<input type="checkbox"/> Vascular/Thin <input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	4	53mm	57mm	56	8	3
Proximate Linear Cutter 75mm	TLC75 TCD75 TCT75	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Regular/Thick <input checked="" type="checkbox"/> Thick	4	73mm	77mm	76	8	3
Proximate Linear Cutter 100mm	TLC10 TCT10	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	4	98mm	102mm	100	8	3



Description	Reload code	Tissue	Closed staple height	Wire diameter	Open staple height	Qty/box
Proximate Linear Cutter 55mm Reload	TVR55 TCR55 TRT55	<input type="checkbox"/> Vascular/Thin <input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	1.0mm 1.5mm 2.0mm	0.20mm 0.20mm 0.23mm	2.6mm 3.8mm 4.5mm	12
Proximate Linear Cutter 75mm Reload	TCR75 TRD75 TRT75	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Regular/Thick <input checked="" type="checkbox"/> Thick	1.5mm 1.8mm 2.0mm	0.20mm 0.20mm 0.23mm	3.8mm 4.2mm 4.5mm	12
Proximate Linear Cutter 100mm Reload	TCR10 TRT10	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	1.5mm 2.0mm	0.20mm 0.23mm	3.8mm 4.5mm	12

Device may be reloaded seven times for a maximum of eight firings.

The third-party trademarks used herein are trademarks of their respective owners.

For complete product information, including full steps for use, indications, contraindications, warnings and precautions, please see Instruction for Use.