

## Absorbable Suture Clip Applier


### Clip the end of sutures in seconds

Compared to tying sutures, the time-saving LAPRA-TY® Absorbable Suture Clip Applier makes clipping Coated VICRYL® (polyglactin 910) Suture as easy as squeezing the ergonomically designed handle.


This reusable device, compatible with LAPRA-TY® Absorbable Suture Clips, anchors and secures each end of a single strand of size 2-0, 3-0 or 4-0 Coated VICRYL Suture for soft tissue approximation for a period of up to 14 days.

### LAPRA-TY Absorbable Suture Clip Applier features

- Clips designed to retain their security on each end of a single strand of Coated VICRYL Suture
- 22° angled jaws designed for maneuverability and visibility for clip placement
- One-piece clips molded from the polymer poly (p-dioxanone)
- 360° shaft rotation
- Can be used in any laparoscopic procedure



KA200



Precision tip securely holds each clip during manipulation and placement.

XC200



Six absorbable clips in each easy-load cartridge, six cartridges per box.

## Ordering Information

Product Code	Description	For Use With	Qty/box
KA200	LAPRA-TY® Absorbable Suture Clip Applier, 33cm length	10mm or 11mm diameter trocars	1
XC200	LAPRA-TY® Absorbable Suture Clips, 6 clips per cartridge, 6 cartridges per box	2-0, 3-0, or 4-0 Coated VICRYL® Suture	6

## How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

### Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

#### **Johnson & Johnson Gateway: [www.jnjgateway.com/commerce](http://www.jnjgateway.com/commerce)**

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

#### **Global Healthcare Exchange (GHX): [www.ghx.com](http://www.ghx.com)**

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

#### **Electronic Data Interchange (EDI): JJHCS EDI Help Line: 1-800-262-2888**

### Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

## Customer Support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at **customersupport@eesus.jnj.com** or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.