

# MICROFLEX®

## COMFORTGRIP®



NATURAL  
LATEX  
NON-STERILE  
AMBIDEXTROUS  
FULLY TEXTURED  
BEADED CUFF  
FOR SINGLE USE ONLY

### POWDER-FREE EXAMINATION GLOVES

Backed by the industry's leading quality system, ComfortGrip® gloves offer the best combination of quality and value in a latex glove. The fully textured surface provides a reliable grip for your daily tasks.

- A great value in a latex glove
- Secure grip
- Comfortable for everyday use
- Fully textured for a great grip



PRODUCT CODE		PACKAGING	
XS	CFG-900-XS	100	Gloves/Box
S	CFG-900-S	10	Boxes/Case
M	CFG-900-M		
L	CFG-900-L		
XL	CFG-900-XL		

GLOVE PROFILE		ASTM D 3578	COMFORTGRIP®
Length	mm (inches)	min. 220 (8.7)	245 (9.6)
Cuff Thickness	mm (mils)	n/a	0.10 (3.9)
Palm Thickness	mm (mils)	min. 0.08 (3.1)	0.13 (5.1)
Finger Thickness	mm (mils)	min. 0.08 (3.1)	0.16 (6.3)
Tensile Strength	Before Aging	min. 18	22
	After Aging	min. 14	20
Elasticity/Elongation	Before Aging	min. 650	700
	After Aging	min. 500	600
Pinhole Defect Rate*	(AQL)	max. 2.5	1.5

\* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED  
INDUSTRIES

Caution: This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.