

MICROFLEX®

▲ EVOLUTION ONE®



NATURAL
LATEX
NON-STERILE
AMBIDEXTROUS
FULLY TEXTURED
BEADED CUFF
FOR SINGLE USE ONLY

POWDER-FREE EXAMINATION GLOVES

Made with premium latex, Evolution One® gloves provide maximum comfort along with enhanced tactile sensitivity. Strong and durable with a secure grip, Evolution One® gloves slip on easily, wet or dry, and provide a consistent performance in mid-range environments.



- Best premium latex offering
- Maximum comfort allows for extended wear
- Slip on easily in wet or dry conditions
- Fully textured for a great grip

OUR BEST PREMIUM LATEX OFFERING

PRODUCT CODE		PACKAGING	
XS	EV-2050-XS	100	Gloves/Box
S	EV-2050-S	10	Boxes/Case
M	EV-2050-M		
L	EV-2050-L		
XL	EV-2050-XL		

GLOVE PROFILE		ASTM D 3578	EVOLUTION ONE®
Length	<i>mm (inches)</i>	min. 220 (8.7)	245 (9.6)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.10 (3.9)
Palm Thickness	<i>mm (mils)</i>	min. 0.08 (3.1)	0.14 (5.5)
Finger Thickness	<i>mm (mils)</i>	min. 0.08 (3.1)	0.15 (5.9)
Tensile Strength	<i>Before Aging</i>	min. 18	24
	<i>After Aging</i>	min. 14	18
Elasticity/Elongation	<i>Before Aging</i>	min. 650	700
	<i>After Aging</i>	min. 500	600
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED INDUSTRIES

Caution: This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.