

MICROFLEX®

# BECAUSE SIZE MATTERS

The first ambidextrous  
exam glove available  
in ½ sizes.

Wearing exam gloves  
that are too loose or  
too tight may lead to  
hand fatigue and possibly  
career impacting injuries.

Your hands are special.  
Your gloves should be too.

**Microflex® Ultraform®** comes in a  
range of ½ sizes for a perfect fit,  
because no two hands are exactly  
the same.

If your glove isn't **Microflex Ultraform**,  
then it isn't half as good.

PLACE YOUR HAND HERE TO  
APPROXIMATE YOUR SIZE.

REQUEST A SAMPLE TO KNOW  
YOUR EXACT SIZE.

Extra-Small 77 mm **NEW!**  
Extra-Small/Small 82 mm  
Small 88 mm  
Small/Medium 92 mm **NEW!**  
Medium 97 mm  
Medium/Large 102 mm **NEW!**  
Large 106 mm  
Extra-Large 115 mm



Microflex®  
**Ultraform®**

**THE BEST FIT** JUST GOT BETTER  
Visit **Microflex.com** for a sample


**Ansell**

## Microflex® Ultraform® Powder-Free Nitrile Examination Gloves

Microflex Ultraform with ERGOFORM technology, an innovation that reduces stress on joints, ligaments and tendons during repetitive tasks. Users experience more comfort and complete tasks with less fatigue minimizing the risk factors that lead to workplace injuries.

- Amazing fit and feel
- Now available in ½ sizes
- Ergonomic Certification
- Bulk Packed 300 gloves/box\*

\*XL 250 gloves/box



## ERGOFORM™ Ergonomic Design Technology



### PRODUCT INFORMATION

Material	Powder-Free Nitrile
Sterility	Non-Sterile
Usage	Single Use Only
Color	Cobalt
Cuff Length	Standard 9.5"
Finish	Textured Finger Tips

### ORDERING INFORMATION

Size	Part #	Packaging
XS	UF-524-XS	300 gloves per box
XS/S	UF-524-XS/S <b>NEW!</b>	
S	UF-524-S	
S/M	UF-524-S/M <b>NEW!</b>	10 boxes per carton
M	UF-524-M	
M/L	UF-524-M/L <b>NEW!</b>	
L	UF-524-L	250 gloves per box
XL	UF-524-XL	

Visit [Microflex.com](http://Microflex.com) to learn more or request FREE samples