



## Introducing: The OmniTrust™ #113 Series Latex Powder Free PL Examination Glove

Our OmniTrust™ Latex Powder Free PL is a reduced protein glove. It has our famous (micro-texture) matte finish for added grip control without loss of tactile sensitivity. This glove is manufactured with a non-chlorination process thus maintaining the stretch characteristics of latex. A unique polymer coating gives this glove excellent donning characteristics. Omni International Latex gloves are flat packed for optimum dispensing in sizes small through extra large.

### Features/Benefits:

- Polymer Coated for Easy Donning
- Micro-Textured for Superior Grip Control
- Non-Chlorinated Natural Rubber Latex
- Ambidextrous, Beaded Cuff, Single Use Only

Order #	Size	Gloves/Box*	Boxes/Case
113-01	S	100	10
113-02	M	100	10
113-03	L	100	10
113-04	XL	100	10

\*Packaged by Weight



**OMNI INTERNATIONAL CORP**  
A Blue Sail Company

360 Route 101, Building 2 • Bedford, NH 03110 • Tel: 888-999-6664 • [info@omnigloves.com](mailto:info@omnigloves.com)