



Introducing: The OmniTrust™ #412 Series Vinyl Powder Free Examination Glove

This glove is made from the highest quality polyvinyl chloride material and formulation available. Our unique Comfort Formulation allows the glove to fit more comfortably to your hand than traditional vinyl gloves. The 412 series is an excellent alternative to latex gloves and all of our gloves are produced with DINP, a safe alternative to the harmful DOP/DEHP plasticizer found in other vinyl gloves. All of this adds up to a superior vinyl glove in quality, comfort, and performance.

Features/Benefits:

- 100% Synthetic Stretch Vinyl (Comfort Formulated)
- Completely Latex FREE (No risk for latex allergy)
- Manufactured without DOP / DEHP
- Smooth Finish for Improved Feel/Tactile Sensitivity
- Ambidextrous Beaded Cuff, Single Use Only

Order #	Size	Gloves/Box*	Boxes/Case
412-01	S	100	10
412-02	M	100	10
412-03	L	100	10
412-04	XL	100	10

*Packaged by Weight



OMNI INTERNATIONAL CORP
A Blue Sail Company