

Compatible Syringes				
The Autoject 2 may be used with the following syringes:				
Syringe Brand	Size	Needle Length	Depth Adjuster	Depth Adjuster Setting (6mm/ 1/4")
B. Braun	1.0ml (1cc)	8mm (5/16")	Long (blue scale)	c
B. Braun	0.5ml (1/2cc)	8mm (5/16")	Long (blue scale)	n
B. Braun	1.0ml (1cc)	12mm (1/2")	Long (blue scale)	e
B. Braun	0.5ml (1/2cc)	12mm (1/2")	Long (blue scale)	q
BD Micro-Fine®/Ultra-Fine®	1.0ml (1cc)	12.7mm (1/2")	Long (blue scale)	c
BD Micro-Fine®/Ultra-Fine®	0.5ml (1/2cc)	8mm (5/16")	Short (black scale)	m
BD Micro-Fine®/Ultra-Fine®	0.5ml (1/2cc)	12.7mm (1/2")	Short (black scale)	o
BD Micro-Fine®/Ultra-Fine®	0.3ml (3/10cc)	12.7mm (1/2")	Short (black scale)	j
Precision Sure-Dose™/Terumo®	1.0ml (1cc)	12mm (1/2")	Long (blue scale)	b
Precision Sure-Dose™/Terumo®	0.5ml (1/2cc)	12mm (1/2")	Short (black scale)	l
Precision Sure-Dose™/Terumo®	0.3ml (3/10cc)	12mm (1/2")	Short (black scale)	l

