



**OWEN MUMFORD**  
Making a World of Difference

**THE AMERICAS**

Owen Mumford Inc.  
1755 West Oak Commons Court  
Marietta, Georgia 30062, USA  
T : (800) 421-6936 F : (770) 977-2866  
E : [info@owenmumfordinc.com](mailto:info@owenmumfordinc.com)

**HEAD OFFICE**

**INTERNATIONAL SALES & MARKETING**

Owen Mumford Ltd.  
Brook Hill, Woodstock  
Oxfordshire OX20 1TU, United Kingdom  
T : 01993 812021 F : 01993 813466  
E : [info@owenmumford.co.uk](mailto:info@owenmumford.co.uk)

**UK / IRELAND SALES & MARKETING**

T : 01608 645555 F : 01608 646546  
E : [info@owenmumford.co.uk](mailto:info@owenmumford.co.uk)

**FRANCE**

T : 02 32 51 88 70 F : 02 32 21 96 95  
E : [contact@owenmumford-fr.com](mailto:contact@owenmumford-fr.com)

**GERMANY**

T : 06026/97750 F : 06026/977529  
E : [info@owen-mumford.de](mailto:info@owen-mumford.de)

**CHINA**

T : 021 6032 3528 F : 021 6032 3782  
E : [info@owenmumford.co.uk](mailto:info@owenmumford.co.uk)

[www.owenmumfordinc.com](http://www.owenmumfordinc.com)

**OUR ACCREDITATIONS INCLUDE:**

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMDCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001




**INNOVATION**  **EXCELLENCE**

UNIFINE/IT/OM/0114/1/US



**A Smarter Pen Needle**

**Unifine<sup>®</sup> Pentips<sup>®</sup>**  
**Pen Needles**



 **OWEN MUMFORD**


# Unifine<sup>®</sup> Pentips<sup>®</sup>

## Exceeding Expectations for Quality and Performance

For over 60 years, Owen Mumford has established a global reputation for innovation and excellence in medical device manufacturing. Produced under strict quality standards, **Unifine<sup>®</sup> Pentips<sup>®</sup>** undergo 100% inspection to ensure our pen needles offer the highest quality and performance.

### Sterile Chamber

Helps keep new needles sterile and protected until ready for use



### Designed for Comfort

Innovative design reduces penetration and glide force for a more comfortable injection

### Universal Compatibility

Fits all major brand insulin and non-insulin pens


## Product Range

 <b>4mm</b> 32G x 4mm	 <b>5mm</b> Mini 31G x 5mm	 <b>6mm</b> Ultra Short 31G x 6mm	 <b>8mm</b> Short 31G x 8mm	 <b>12mm</b> Original 29G x 12mm
100ct 08470-3540-01 30ct 08470-1140-01	100ct 08470-3550-01 30ct 08470-1150-01	100ct 08470-3590-01 30ct 08470-1190-01	100ct 08470-3530-01 30ct 08470-1130-01	100ct 08470-3529-01


## The Unifine<sup>®</sup> Pentips<sup>®</sup> Advantage

- Triple sharpened for lower penetration force
- Anti-coring technology for improved drug delivery
- Advanced silicone for a lower glide force
- Thin wall technology for improved drug flow and low delivery force

## Proven Performance



Unifine<sup>®</sup> Pentips<sup>®</sup> With Anti-Coring



Other Pen Needle Without Anti-Coring

## A Smarter Pen Needle

Covered by most third party payors (including Medicaid and Medicare), Unifine<sup>®</sup> Pentips<sup>®</sup> offer a premium pen needle without the premium price. Patients will love the comfortable injection experience that comes with lower glide and penetration force.

Pharmacists appreciate how competitive price and generous reimbursement make Unifine<sup>®</sup> Pentips<sup>®</sup> a smarter business choice over other leading national brands.

