

Protection against CLABSIs*

WHERE YOU
NEED IT,
WHEN YOU
NEED IT.



New Strip Format

Ready access from the IV pole

+Prevantics®

3.15% Chlorhexidine Gluconate (CHG) and
70% Isopropyl Alcohol (IPA)

Device

SWAB

- Disinfects needleless access sites prior to use
- Designed for convenience
- Strip contains 8 swabs

*CLABSI= Central line-associated bloodstream infection



+PDI®
BE THE DIFFERENCE

PATIENT CARE

ENVIRONMENT OF CARE


INTERVENTIONAL CARE

Product features and benefits:

- The first and only FDA-cleared medical device featuring a 3.15% (w/v) Chlorhexidine Gluconate and 70% (v/v) Isopropyl Alcohol formulation
- Confirmed efficacy at the 5 second scrub time and 5 second dry time
- Achieved a $>4.0 \log_{10}$ reduction versus a leading alcohol impregnated cap¹ against 7 tested clinically relevant microorganisms known to cause CLABSIs²



Intuitive, easy-to-use prep-pad format.



Meets evidence-based clinical guidelines and recommendations from the following organizations:



	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Prevantics® Device						
Swab	B19600	n/a	10/100's	4.30 lbs	0.28 ft	21/7
Swab Strip	B123ST	n/a	10/160's	4.68 lbs	0.518 ft	8/5

For a free sample, visit pdihc.com/PrevanticsDeviceStrip

References: 1. Data on file, study report 140304-250. 2. CDC HAI Prevalence Survey. Magill SS, Edwards JR, Bamberg W, et al. Multistate Point-Prevalence Survey of Health Care–Associated Infections. *N Engl J Med.* 2014;370:1198-208.