



SANI-CLOTH® BLEACH

GERMICIDAL DISPOSABLE WIPE

Ideal for disinfecting high risk areas endemic with Multidrug-Resistant Organisms, *Clostridium difficile* spores and Norovirus.

Sani-Cloth® Bleach Benefits:

- Meets CDC, OSHA and CMS Tag F441 guidelines
- Bactericidal, Fungicidal, Tuberculocidal, Virucidal
- Compatible with a broad range of surfaces and equipment in healthcare
- For use when taking a pathogenic specific approach to disinfection

EPA Reg. No. 9480-8

Contact Time
(Minutes)



of Microorganisms
Effective Against

50

Disinfection
Formulation

**1:10
Bleach
Dilution**

Effective against 50 microorganisms in 4 minutes including the following MDROs, bloodborne pathogens and viruses:¹

- *C. difficile* spores
- *Acinetobacter baumannii*
- *Klebsiella pneumoniae*
- *Candida albicans*
- ESBL-resistant *E. coli*
- ESBL-resistant *Klebsiella pneumoniae*
- Norovirus
- MRSA
- VRE
- VRSA
- HIV
- HBV
- HCV

SANI-CLOTH®

The #1 brand of disinfectant wipes in healthcare.²

¹ See Technical Data Bulletin for a complete list of microorganisms. ² GHX Intelligence Data 2013.



ACCESSIBILITY: Deep well lid provides ample space to store and access next wipe.



COMPLIANCE: Prominent contact time for easy identification and improved compliance with accreditation standards.



IDENTIFICATION: Categorization of pathogens allows for quick identification of relevant efficacy claims.

ORDER INFORMATION

SANI-CLOTH® BLEACH (1:10 BLEACH DILUTION)	Reorder No.	Wipe Size	Packaging	Case Wgt. (lbs.)	Case Cube (ft.)	TI/Hi
A. Large Canister	P54072	6" x 10.5"	75/can, 12 cans/cs	22.65	1.540	10/3
B. Extra-Large Canister	P25784	7.5" x 15"	65/can, 6 cans/cs	14.30	1.012	10/4
C. Large Individual Packets	H58195	5" x 7"	40/box, 10 boxes/cs	8.10	0.505	19/5
D. Extra-Large Individual Packets	U26595	11.5" x 11.75"	40/box, 3 boxes/cs	8.20	0.4998	24/3
E. Pail	P7007P	7.5" x 15"	160/pail, 2 pails/cs	16.5	0.953	10/5
F. Refill	P700RF	7.5" x 15"	160/refill, 2 refills/cs	15.2	0.529	18/3



SANI-CLOTH® EDUCATIONAL & COMPLIANCE ACCESSORIES

Developed to assist you in facility protocol compliance, as well as meet The Joint Commission and CMS requirements.^{3,4}



PDI® 3-in-1 Universal Sani-Bracket®
for wall/pole mounting of a
Sani-Cloth® Canister



Sani-Canister Caddy®
for placement on flat
surfaces (bilingual)



Sani-TAG™ Equipment ID System
for identification of
clean/soiled equipment



Educational wall charts and signage
(bilingual)

www.pdihc.com
1-800-999-6423

© 2014 PDI. All rights reserved.
PDI01147161

³Infection Prevention and Control (IC) Standard IC.02.02.01 and Elements of Performance (EP) 1 and 2, The Joint Commission, October 2009. ⁴Guideline for Disinfection and Sterilization in Healthcare Facilities, Centers for Disease Control and Prevention, 2008.

PDI
BE THE DIFFERENCE™