



Antimicrobial keypad, handrails
and transport handle

FALL PREVENTION

Safe patient handling is a primary concern for healthcare providers

- ✔ **“Live” Handrails** allow the patient to hold onto handrails for stability while the scale accurately calculates weight
- ✔ **Low Profile Platform** requires only minimal step up, offering the patient a more comfortable and secure weighing experience
- ✔ **Heavy-Duty Construction** and low center of gravity comforts unsteady patients and reduces the chance of tipping over

INFECTION PREVENTION

Reducing Healthcare Associated Infections is a major challenge for medical facilities

- ✔ **Antimicrobial*** keypad, handrails and transport handle help reduce patient and staff exposure to various pathogens
- ✔ **Contemporary Design** eliminates places for dirt and contaminants to collect

DURABILITY & EASE OF USE

Facilities require products that withstand frequent use, movement, and disinfection

- ✔ **Convenient Transport Handle** allows users to easily move the scale throughout the facility
- ✔ **Impact Resistant and Floor Protective Wheels** ensure smooth movement and enhanced portability
- ✔ **Simple User Interface** eliminates confusion so caregivers can focus on getting accurate weight while keeping patients safe

*Antimicrobial protection is limited to the treated article and does not protect a user against disease causing bacteria



“Live” handrails for patient stability and accurate weight measurement



SPECIFICATIONS



- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 19-1/4" (w) x 21-3/4" (d) x 1-3/4" (h) / 489 mm (w) x 552 mm (d) x 44 mm (h)
- Handrail Width: 21-1/2" / 546mm
- Display: 3/4" LCD
- 6" x 1-1/2" Rubber Wheels: 2
- Power Source: 6 D batteries and 110-240V adapter included
- Connectivity: USB
- Functions: LB / KG Conversion, LB / KG Lock Out, Hold, Recall, Reweigh, Auto Zero, Auto Off, Audible Signal, Low Battery Indicator
- 2 Year Limited Warranty



DESIGNED AND
MANUFACTURED
IN THE . . .

USA