



9990 NW 14<sup>th</sup> Street,  
Unit #105  
Doral, Florida 33172  
Tel: 305.503.2729  
Fax: 786.664.3347  
info@ezerusa.com  
www.ezerusa.com

Distributed by



# ECP-5200

## CHART PROJECTOR



### Accessories



Polarized Glasses

Screen

### Optional



Projector Stand

Projector Wall Mount

# ECP-5200

## SPECIFICATION

Projection Distance:	1.5 - 10m
Project Magnification:	30° x (In 5m refraction)
Number of Charts:	33
Mask & Filter:	1 Open, 5 Horizontal lines, 5 Vertical lines, 15 Single isolation, 1 Red/Green
Power Supply:	100 - 240V AC 50/60Hz
Lamp:	LED Lamp
Power Consumption:	60W
Power Saving:	5 min
Tilt Angle:	+/- 20°
Dimension & Weight:	197 (W) x 335 (D) x 163 (H)mm, 5.6kg
Accessories:	Remote Control / Power Cable / Dust Cover / Battery / Fuse/ Spare Lamp/ Table Stand / Screen
Optional Accessories:	Wall Mount Bracket / RS-232C (Interface) / Unit Adaptor / Red-Green Filter / Polarized Filter

Designs and details can be changed without prior notice for the purpose of improvement.



# ECP-5200

## CHART PROJECTOR



### Digital Practice

ECP-5200 can be easily connected to the EDR-5200 (Digital Refractor) for fast and accurate results, minimizing transcript errors.

### Digital Practice



### Rapid and Exact Chart change

This feature provides you easy and accurate measurement.

### Compact Design

With a compact and light design, the ECP-5200 can be placed in a small and narrow space.

### Wide Range of Measurement Functions

Various measurement modes include: binocular balancing, strabismus, aniseikonia, stereovision, binocular fusion, suppression, etc. With multiple charts: Red/Green polarizing filter and Maddox Rod to check heterophoria, we are able to achieve higher reliability of vision tests.

### Clear Image

ECP-5200: long lasting white projection light. Provides an effective and accurate measuring environment, with a sharp, clear projection in high resolution.

E 400	S 400	200	3ME 200	ME 100
FRA 300	SB2 300	AW 100	MBE 100	MBE 100
HDN 150	386 150	70	MBE 100	MBE 100
CVYHO 100	8888 100	W-A 30	MBE 100	MBE 100
LFPAZ 30	8888 80	S-V 30	MBE 100	MBE 100
OTCN 30	8888 70	Ma-V 20	MBE 100	MBE 100
TMFOZ 50	2888 60		MBE 100	MBE 100
NETOP 30	8888 50		MBE 100	MBE 100
ALVAH 40	8888 40		MBE 100	MBE 100
HVNAO 30	8888 30		MBE 100	MBE 100
LZCR 15	8888 20		MBE 100	MBE 100
CTFDK 20	8888 20		MBE 100	MBE 100
YOVAN 15	8888 15		MBE 100	MBE 100
LOHNT 15	2888 15		MBE 100	MBE 100
KEZFR 10	8888 10		MBE 100	MBE 100
			MBE 100	MBE 100
			MBE 100	MBE 100

20/20 Chart

### Wide Range of Charts

The ECP-5200 with 23 charts including Red/Green and Polarizing filters are the new solution for Binocular balance test, Stereo test, Aniseikonia test, Fusion test.



### Various Control Methods

Charts can be changed from the Ezer EDR-5200 (Digital Refractor) or with the remote control, conveniently from your desk.