

DRE WS350

Diagnostic Wall Mount System

Equipment for the way *you* operate

The DRE WS350 Wall Mount Diagnostic System includes all the instruments you'll need to perform a complete diagnostic examination: illuminating scopes, thermometers, specula and more. The system is designed for accuracy and ease of use.

Key Components:

- **Wall Mount Transformer with Handles**
 - » Comes Standard with Two Diagnostic Handles
 - » Spiral Cords Expandable up to 3 meters
- **3.5V Coaxial Ophthalmoscope**
 - » Enables easy access to undilated pupils for accurate diagnosis
 - » Halogen lamp provides light output for true tissue color and consistent, long-lasting illumination
- **3.5V Fiber Optic Otoscope**
 - » >2X Magnification
 - » Swiveling wide-angle viewing lens allows for instrumentation during magnification
 - » Sealed system for pneumatic otoscopy and accurate diagnosis of otitis media
- **Ear Thermometer**
 - » Preheated probe tip and sensor
- **Sphygmomanometer**
 - » Dual scale dial face provides unmatched accuracy
 - » Precision mechanism mounted on a special support to ensure complete stability of the measuring system and the zero-point
 - » Latex and Mercury Free
- **Reusable Ear Specula**
 - » 2.7mm & 4.2mm, 35 pieces each
 - » Smooth finish
 - » Autoclavable
 - » Positive locking design
- **Specula Dispenser-Tip Dispenser for all Specifications of Ear Specula**
- **Wooden Plate mount (large size)**

