

# DRE Velocity EKG

12-channel PC-based EKG

Equipment for the way *you* operate

## Portable, affordable PC-based EKG solution with powerful data and analysis functions

PC, printer and stand not included



### Affordable EKG Solution

The DRE Velocity EKG PC-based EKG includes a data acquisition module and analysis software for acquiring, analyzing, storing and reviewing 12-channel ECGs. It's great for hospitals, physician offices and laboratories.

### Easy to install!

Simply install the software and connect the DRE Velocity EKG to your PC or laptop and you're ready to go.

### Features include:

- » 12-lead real-time waveforms and ST segment/trend graphics display.
- » The software's powerful data and analysis functions assist with diagnosis
- » Email reports in XML, SCP, PDF, Word and JPEG formats.
- » Advanced analysis functions: Vector ECG (3-D), timed VCG, signal average electrocardiograph (SAECG), heart rate turbulence (HRT), QT interval dispersion, high frequency ECG, and frequency spectrum cardiogram.
- » Optional customized exercise protocol for Bruce and Modified Bruce control treadmill; bicycle and exercise blood pressure devices compatible with leading treadmills.
- » Print clear reports the moment you need them on plain paper.



HF ECG

FEKG

Powerful database

VCG 3-D



SAECG

HRV

HRT

QT interval dispersion

### Options

- » Stress test and exercise blood pressure monitoring function
- » Support manual input of blood pressure value



# DRE Velocity EKG

12-channel PC-based EKG

Equipment for the way *you* operate

## Technical specifications

Dimensions	5.8" (l) x 3.9" (w) x 1.6" (h) 148 mm (l) x 100 mm (w) x 40 mm (h)
Weight	Approx. 0.46" (210 g)
Connection to PC	USB port
Power supply	DC 5V from USB port
Input circuit	Defibrillation protection built-in
Lead	Standard 12 leads
Acquisition mode	12 leads simultaneously
Sampling rate	1000Hz
A/D converter	12 bits
Input voltage range	$\pm 5\text{mVpp}$
Time constant	$> 3.2\text{s}$
Frequency response	0.05Hz ~ 150 Hz
Calibration voltage	1mV $\pm$ 2%
Sensitivity	2.5, 5, 10, 20 (mm/mV)
Input impedance	50M $\Omega$ (10Hz)
Input circuit current	$< 50\text{nA}$
Noise level	$< 15 \mu\text{Vp-p}$
Baseline drift filter	Baseline forward and backward filter (patent pending)
Patient leakage current	$< 10 \mu\text{A}$ (220V - 240V)
Patient auxiliary current	$< 0.1 \mu\text{A}$ (DC)
Filter:	
- EMG filter	25Hz/35Hz/45Hz/off
- Baseline drift filter (patent pending):	Strong/weak/off
- Lowpass filter	150 Hz/100 Hz/75 Hz
- AC filter	50 Hz/ 60 Hz
CMRR	$> 110\text{dB}$
Rhythm lead	Single lead or three leads selectable
Safety standard	IEC II/CG

## PC Requirements

Operating system	Windows 2000, Windows XP, Windows Vista or Windows 7
CPU	Pentium P4, Celeron D 310 or above
RAM	256 MB or above
HDD	40G or above
Interface	At least 2 USB ports free

