

DRE Millennium 5

Hospital Stretcher

Available with easy-to-steer fifth wheel system

Benefits

- ✓ High-low adjustable mattress range is 13.5"
- ✓ Back-rest lifting range is 0 - 75° and controlled by oil hydraulics
- ✓ Knee-rest lifting range is 0 - 30° and controlled by adjustable crank
- ✓ Central braking system with steering function
- ✓ Trendelenburg range is $\pm 18^\circ$
- ✓ Two IV poles (one at the head and foot)
- ✓ 4" high density PU foam mattress covered with synthetic leather
- ✓ Collapsible side rails on both sides that can fold manually
- ✓ 8" steering wheel central brake on four casters

Features

- ✓ 4" x 28.5" x 71.5" Pad
- ✓ 8" Casters
- ✓ 27" Wide Patient Surface
- ✓ Fold-Down Side Rails; Side Rails 14.5" x 39"
- ✓ Stationary Foot Section
- ✓ Four Brake/Steer Pedals
- ✓ Dual Pedestal Hydraulics
- ✓ Dual Side-Mounted Foot Control Hydraulics
- ✓ Hands-Free Trendelenburg and Reverse Trendelenburg
- ✓ Steel bed frame with powder coating
- ✓ Mattress base made of ABS plastic board
- ✓ Bumpers made of durable ABS plastic and located at each corner
- ✓ Base cover made of ABS plastic with space for an oxygen bottle and emergency devices



SPECIFICATIONS



Weight Capacity: 550 lbs



Overall Length: 83"

Overall Width: 33"

Mattress Surface:
71.5" (L) x 28.5" (W)

Height Adjustment Range:
22.5" - 36.5" with 8" Casters



**Trendelenburg/
Reverse Trendelenburg:** $\pm 18^\circ$

Backrest Adjustments: +90°



Casters:
8" with central locking one tracking 4 brake/steer pedals

I.V. Attachment Sockets:
2

The fifth wheel system saves power and allows for easy maneuverability

397-01-2017-04