

# DRE Lucerne 360

## Universal Operating Table


### Effective, economical choice for your surgery suite

With the DRE Lucerne 360, you get premium features for as much as 30% less than the cost of competing surgery tables. Crafted for quality, the Lucerne 360 provides intuitive movement and stability, setting new standards in patient positioning. It features a C-Arm compatible table top with a 225 kg (500 lb) lift capacity.

The Lucerne 360 has many advanced positioning options, such as kidney bridge positioning, 360° horizontal rotation (180° each direction) and a horizontal slide motion that helps create a greater field of view for C-arm use.

#### Features

- ✓ C-arm interface capabilities to meet imaging requirements
- ✓ Silent electro-hydraulic movements
- ✓ Mobile anti-static castors with a floor lock
- ✓ Remote control provides high operation safety and precise patient positioning with a safety lock
- ✓ Easy access to electro-hydraulic system for service
- ✓ Kidney bridge and flex positioning



345-01-2017-04

#### SPECIFICATIONS



**Maximum lift capacity:**  
500 lb. (225kg)



**Bed size:**  
80.7" x 19.6"  
(2,050 mm x 500 mm)

**Height:**  
26.3" - 43.3"  
(670 mm x 1,100 mm)



**Horizontal Rotation:**  
Full 360°

**Vertical Rotation:**  
0 ~ 180° each direction

**Reverse Trendelenburg:**  
Antersion: 30°  
Retroversion: ≥30°

**Back Plate Adjustment:**  
+75° ~ -20°

**Head Plate Adjustment:**  
40°/90°

**Leg Plate Adjustment:**  
90°

**Leg Decline:**  
105°



**Voltage:** 220V ± 22V

**Frequency:** 50Hz ± 1Hz