

# DRE Torino 550

## Mobile Surgery Table

### Multi-purpose positioning capabilities

The DRE Torino 550 is a multi-purpose mobile surgery table with electronic positioning capabilities for a variety of surgical procedures. Using the handheld remote control, operators can adjust the DRE Torino 550 into a variety of positions including Trendelenburg and reverse Trendelenburg. Lateral tilt and longitudinal sliding are also standard features. The Torino 550 has large, durable casters that allow clinicians to move patients seamlessly through pre-op, to the surgical suite and into post-operative care.

#### Features

- ✓ Longitudinal slide up to 350 mm for easy C-arm access.
- ✓ Table weight load capacity rated at 550 pounds.
- ✓ Adjustable height from 28 to 40 inches for good access to the operative field.
- ✓ Table movements include: Trendelenburg, reverse Trendelenburg, leg and back plate adjustment and lateral tilt.
- ✓ Large casters and convenient braking system for easy and convenient maneuvering throughout any surgical facility.
- ✓ Battery backup for continued use in the event of a power loss.
- ✓ Manual leg rest.
- ✓ Corded remote control panel for exact positioning.
- ✓ Comfortable memory padding is waterproof and antistatic.
- ✓ Added safety with override panel controls on the base of the unit.
- ✓ Optional — Kidney bridge, shoulder supports, arm supports, long fixing clamps, leg and foot supports and more.



509-2017-11-10

#### SPECIFICATIONS



**Full Length:**  
81.1in (2060mm)

**Full Width:**  
21.7in (550mm)

**Tabletop Elevation:**  
28in (710mm) - 40in (1010m)

**Kidney Bridge Elevation:**  
4.3in (110mm)

**Horizontal Sliding:**  
13.8in (350mm)



**Head Section:**  
45° up - 90° down

**Back Section:**  
80° up - 20° down

**Lateral Tilt:**  
15° right, 15° left

**Trendelenburg:**  
25°

**Reverse Trendelenburg:**  
25°