

DRE Versailles P100

Mobile Surgery Table

The complete multi-function mobile operating room table

The DRE Versailles P100 is more than a surgery table. It lets you move patients from pre-op to the operating room to post-op without transferring the patient from a stretcher to a surgery table and back. The DRE Versailles P100 does it all by incorporating the mobility of a stretcher into a fully-functional surgery table.

The powered DRE Versailles P100 operating room table is more mobile than most stretchers or surgical tables and is a total solution for hospitals and surgery centers. This multi-function table can replace several devices at once to create an extreme value. Rely on the DRE Versailles P100 for exceptional stability, full C-arm access, patient transport, and complete surgical positioning.

Features

- ✓ Traversing patient platform
- ✓ Full C-arm coverage
- ✓ Superior surgical access
- ✓ Lateral tilt
- ✓ Excellent stability
- ✓ 550-lb. (250 kg) weight limit
- ✓ Four remote-activated powered functions



593-2017-11-21

SPECIFICATIONS



Max. Patient Weight Limit:
550 lbs (250 kg)

Trolley Weight:
352 lbs (160 kg)

Leg Section Weight:
13.2 lbs (6 kg)



Height:
23.8" - 39.5"
(605 - 1,005 mm)

Trolley Length:
82.7" (2,100 mm)

Side Rail Width:
25.8" (655 mm)

Cotside Width:
30.3" (770 mm)

Brake Width (Brakes off):
33" (835 mm)

Mattress Dimensions:
80.3" x 23.6" x 3"
(2,040 x 600 x 65mm)



Trendelenburg/Reverse Trendelenburg:
20°/12°

Lateral Tilt:
+12°/-12°

Backrest Articulation:
0 - 65°

Head Section Articulation:
+30°/-30°

Transport and Recovery

Patient safety rails fitted for patient comfort and safety

Unique retractable steering mechanism engaged for maximum maneuverability

Ophthalmology

Patient platform traversed head end for access

Ophthalmic head section fitted



Endoscopy

Designed for every stage of endoscopy — carries patients from pre-op to procedure to recovery.

Offers C-arm compatibility for imaging and the benefits of a fully-featured surgical table.

C-arm access

Traversing top allows for full C-arm access



Gynecology Features

Patient platform traversed to foot end for access

Leg section removed for access

Head section removed for anesthesia access

Fitted with optional powered leg supports



Shoulder Arthroscopy Features

Patient platform traversed foot end for access

Optional arthroscopy backrest fitted



Additional Trolley Features

- ✓ Powered Functions
- ✓ Lateral Tilt
- ✓ Head End Traverse
- ✓ Foot End Traverse
- ✓ K8 Pressure Care Mattress
- ✓ Removable Ultra Lightweight Leg Section
- ✓ Removable and Rotating Patient Safety Rails

- ✓ Push Button Release and Self Locking Controls
- ✓ 5th Wheel Steering Mechanism
- ✓ Advanced 4 Caster Braking System
- ✓ E Size Oxygen Cylinder Trough (Side Loading)
- ✓ Transfusion Pole with Mounting Position
- ✓ 150 mm Caster