

Multi-Use, Ultrasound Table with Stirrups

- Three-piece top with pneumatic adjustable backrest, manual footrest, and stirrups
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails, locking casters & more
(page 20)

80069	length	width	height range
overall	76"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
	193 cm	68.58 cm	54.61 cm – 95.25 cm
backrest	38"	27"	
	96.52 cm	68.58 cm	
footrest	12"	18.5"	
	30.48 cm	47 cm	



80069



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters