

DRE SureVac 22L

Vacuum Autoclave

Dependable sterilization at an affordable price

The DRE SureVac 22L Vacuum Autoclave is a 22-liter unit ideal for busy facilities in need of simple, economical sterilization. With an intuitive user interface and automatic self-diagnostic capabilities, the DRE SureVac 22L provides safe and efficient sterilization every time.

Features

- ✓ Easy-to-use, 22-liter vacuum autoclave.
- ✓ Microprocessor controlled unit features an intuitive user interface for safe and efficient operation.
- ✓ Digital display indicates parameters and conditions.
- ✓ Automatic self-checking mechanism diagnoses and protects against overheating or over-pressurization.
- ✓ Autoclave utilizes cool, filtered air at the end of each cycle to decrease instrument turnaround time.
- ✓ Self-contained operation prevents external steam release during the sterilization process.
- ✓ Interlocking door automatically engages at the beginning of each cycle for a more secure sterilization environment.
- ✓ Fully automatic, controlled quick cycles.
- ✓ Adjustable program parameters.
- ✓ Chamber can fit three to five instrument trays.
- ✓ Comes standard with drain hose and adaptors, filters, autoclave rack and shelves.
- ✓ Available with a thermal printer as an additional option.



496-01-2016-11

SPECIFICATIONS



Autoclave inner chamber:
10 in x 17 in

Outer dimensions:
24 in D x 17.5 in W x 15.5 in H



Frequency: 60Hz

Voltage:
one-phase 110V±11V



Temperature of environment:
41~104°F



Relative Humidity: ≤80%



Atmospheric Pressure:
41~104°F

Sterilization Processes

250°F/universal: Sterilizing time: 20min | Sterilizing temperature: 250°F | Sterilization pressure: 16.0psi | Drying time: 15min

273°F/loose instruments: Sterilizing time: 4 min | Sterilizing temperature: 273°F | Sterilization pressure: 30.5psi | Drying time: 4 min

273°F/universal: Sterilizing time: 4min | Sterilizing temperature: 273°F | Sterilization pressure: 30.5psi | Drying time: 15min

273°F/cotton: Sterilizing time: 4min | Sterilizing temperature: 273°F | Sterilization pressure: 30.5psi | Drying time: 20min

B-D TEST: Sterilizing time: 3.5min | Sterilizing temperature: 273°F | Sterilization pressure: 30.5psi | Drying time: 10min

VACUUM TEST: 15min