



Hy21[®] Bedpan support

Tender description

Hy21 yellow bedpan support designed for Hygienic Covers
Bedpan washer compatible
Recyclable
Made in Canada 

Product code

Case code: BA-BASS-HY21-JAU

Presentation

Hy21[®] yellow bedpan support in polypropylene. It is used in association with an hygienic cover[®] – bedpan bag. It allows urine and stools collection. Institutions usually use the yellow colour to indicate a contaminated patient.

Colour

Yellow.

Support dimensions

Length = 15.5" – (38.5cm). Width = 11.75" - (29.6cm). Height = 4" – (10.6cm).

Composition

Polypropylene. No latex and no phthalates.

Shelf life

No expiry date.

Efficiency

Weight resistance = 300 lbs (136 kg).

Manufacturing process

Injection moulding.

Use

For any person, men or women, who needs a bedpan whether in a health care institution or at home. Particularly recommended for isolation case.

Instructions for use

Refer to the tutorial for more detail.

Precautions for use

Single-patient use. It is recommended to discard the bedpan support after the patient's departure.

Cleaning / Sterilization

The bedpan support must be regularly cleaned and when soiled. It is recommended to use a disinfectant detergent or an enzyme product associated with a mechanical rubbing action. Follow the instructions on disinfectants to maximize their efficiency (dilution, contact time ...). It can also be cleaned in a bedpan washer without exceeding the temperature of 203°F or 95°C.

Disposal

Dispose of product in regular waste or according to the protocol of the institution. Don't dispose the bag or pad in toilet. The support is recyclable (Plastic #5 - PP).

Benefits

Increases hygiene, saves time, reduces odour, simple and risk-free disposal. The kit is a patented product with hole-shaped handles. Limits the spread of pathogens.
 The Hy21[®] bedpan supports can be stacked onto each other to facilitate storage and thereby save space.

Identification

Each case of bedpan supports is identified with a lot number.

Regulatory

FDA: class I medical device exempt from 501(k).

CE marking: Class I medical device according to the rules of annex IX of European Council Directive 93/42/EEC.

Risks / Safety

This product was evaluated and is not considered hazardous based on the hazard criteria established by the OSHA (Occupational Safety and Health Administration). Safety data sheet is available upon request.

Handling/Storage

Store in a clean and dry place, away from sun.

Packaging

Unit	Not packaged individually	38.5 x 29.6 x 10.6 cm 15.5 x 11.75 x 4 in
Case	20 stacked bedpan supports per case	40.5 x 60 x 32 cm / 4.5 kg 16 x 12.5 x 23.75 in / 10 lbs
Pallet	18 cases per pallet or 360 bedpan support	122 x 102 x 133 cm / 104 kg 48 x 40 x 52¼ in / 230 lbs