

GO GLOV® - Disposable Toilet Glove *Pratical and Effective*



- ✓ Rectangular Glove - 16 x 22 cm
- ✓ Weight 75g.
- ✓ 2 sides
- ✓ Very soft fabric, lint-free
- ✓ Fibers meet the requirements of the OE-KOTEX® Standard 100 Class 1
- ✓ Very resistant
- ✓ Can be used on baby's skin
- ✓ 100% recyclable

● Ready to use, the GO GLOV® is used for personal hygiene to patients in medical facilities, nursing homes, home care and maternity. Made in very soft and lint-free cloth, they are non-allergic and 100% recyclable. The fibers meet the requirements of the OE-KOTEX® Standard 100 Class 1, which allows to be use on the skin of babies.

- Stitch bonded fabric : 50% rayon, 50% polyester
- Weight in grams : 75g/m²
- Low deformation:
 - Breakage machine way: 1250 N/m :
 - Breakage machine width : 1250 N/m
 - Absorption capacity: 750%

- Folded fabric and fused with ultrasonic points
- Non allergenic
- Resistance to sterilization:
 - Autoclave : yes
 - Ethylene oxide : yes

To order:

Product code: GA-SEC-01
Case: 1000 gloves

www.hygie.com

Toll Free : 1 (866) 588-2221

Fax : (450) 444 6222

Canada, USA