



Hy21[®] Vomit support

<u>Tender description</u>	Hy21Vomit support designed for Hygienic Covers Bedpan washer compatible Recyclable Made in Canada 
Product code	Case code: HY-SUPP-VOMO-000
Presentation	Polypropylene Hy21 [®] vomit support. It is used in association with an hygienic cover [®] – vomit bag containing a super absorbent pad. It allows vomit collection.
Colour	White.
Support dimensions	Length: 7.5" - (19cm) x Width: 4.5" - (11.5cm) x Height: 5.1" - (13cm).
Composition	Polypropylene. No latex and no phthalates.
Shelf life	No expiry date.
Manufacturing process	Injection moulding
Use	For any person, men or women, experiencing vomiting, whether in a health care institution, at home or during a trip. Particularly recommended for oncology services, dialysis, post-operative recovery room, as well as for emergency services, such as ambulances, etc.
Instructions for use	Refer to the tutorial for more detail.
Precautions for use	Single-patient use. It is recommended to dispose of vomit support after the patient's departure.
Cleaning / Sterilization	The vomit support must be regularly cleaned and when soiled. It is recommended to use a disinfectant detergent or an enzyme product associated with a mechanical rubbing action. Follow the instructions on disinfectants to maximize their efficiency (dilution, contact time ...). It can also be cleaned in a bedpan washer without exceeding the temperature of 203°F or 95°C.
Disposal	Dispose of product in regular waste or according to the protocol of the institution. The support is recyclable (Plastic #5 - PP).
Benefits	Increases hygiene, saves time, reduces odours, simple and risk-free disposal.
Identification	Each case of vomit supports is identified with a lot number.
Regulatory	FDA: class I medical device exempt from 501(k). CE marking: Class I medical device according to the rules of annex IX of European Council Directive 93/42/EEC.
Risks / Safety	This product was evaluated and is not considered hazardous based on the hazard criteria established by the OSHA (Occupational Safety and Health Administration). Safety data sheet is available upon request.
Handling/Storage	Store in a clean and dry place, away from sun.

Packaging	Unit	Not packaged individually	19 X 11.5 X 13 cm / 75.4 g 7.5 x 4.5 x 5.1 in / 0.17 lbs
	Case	20 vomit supports per case	27 x 60 x 22 cm / 2.3 kg 10,5 x 23,5 x 8,5 in / 5 lbs
	Pallet	40 cases or 800 vomit supports	122 x 107 x 120 cm / 113 kg 48 x 42 x 47¼ in / 250 lbs