

### Wall support system for Hygie products

<b>Presentation</b>	Wall support system for boxes of Hygie Hygienic covers® and Hy21® supports. Adjustable to your needs.
<b>Color</b>	White.
<b>Composition</b>	Cold roll steel with a powder coating finish. No latex and no phthalates. No porous material.
<b>Shelf life</b>	No expiry date.
<b>Use</b>	For any healthcare institutions user of Hygie products.
<b>Instructions for use</b>	Hang the support system on the wall. Insert the Hygie box or the Hy21 support in the wall support system. Remove the box's lid along the pre-cut perforations.
<b>Precautions for use</b>	None.
<b>Cleaning / Sterilization</b>	Clean with a disinfectant wipe or with water and soap. Do not submerge in water.
<b>Disposal</b>	Dispose of product in recycling or according to the protocol of the institution.
<b>Benefits</b>	Easy Hygienic Covers dispensing. Hy21 supports always within arm's reach, ready for use.
<b>Identification</b>	Each support is identified with a lot number.
<b>Handling/Storage</b>	Store in a clean and dry place, away from sun.
<b>Models available</b>	4 models



**Product code : HY-HOLD-PART-001**

Can hold 1 box of Bedpan Hygienic Covers OR 1 x 35 unit GOVOM display

Dimensions: 9.75 x 8.75 x 5.25 in - 25 x 22 x 13.3 cm

Case dimensions : 19 x 9.5 x 6 po – 48.2 x 23 x 15.2 cm

Case weight : 1.3 lbs – 0.6 kg



**Product code : HY-HOLD-PART-002**

Can hold 1 box of Bedpan Hygienic Covers AND 1 Hy21 Bedpan support

Dimensions: 14.25 x 9.75 x 5.25 in - 36 x 25 x 13.3 cm

Case dimensions : 19 x 9.5 x 6 po – 48.2 x 23 x 15.2 cm

Case weight : 1.65 lbs – 0.75 kg



**Product code : HY-HOLD-PART-003**

Can hold 1 box of Urinal Hygienic Covers AND 1 box of Vomit Hygienic Cover AND 1 Hy21 Urinal support AND 1 Hy21 Vomit support

Dimensions: 13.25 x 10.75 x 5.5 in - 33.5 x 27.3 x 14 cm

Case dimensions : 19 x 9.5 x 6 po – 48.2 x 23 x 15.2 cm

Case weight : 2.1 lbs – 0.95 kg



**Product code : HY-HOLD-PART-004**

Can hold all the Hygie products

Dimensions: 26.5 x 13.25 x 5.5 in - 67 x 33.5 x 14 cm

Case dimensions : 19 x 9.5 x 6 po – 48.2 x 23 x 15.2 cm

Case weight : 3.2 lbs – 1.45 kg