

The FastPack® IP Free T4 Immunoassay



FAST / SIMPLE / ACCURATE

The FastPack® IP Free T4 Immunoassay

The FastPack® IP Free T4 Immunoassay, for use on the FastPack® IP System, allows any testing facility to produce point-of-care results in 7 minutes with less than one minute of hands-on time.

- Enhanced quality of service for patients
- Ability to monitor and treat in one visit
- Improved patient compliance and outcome
- Laboratory accuracy and precision



SEKISUI
DIAGNOSTICS
Because every result matters™

The FastPack® IP Free T4 Immunoassay

Simple Process



Dispense sample into the FastPack® Injection Port



Place FastPack® into the FastPack® IP Instrument



Press the START button

The FastPack® IP Free T4 Assay Specifications:

- Quantitative measurement of Free T4
- Lower Limit of 0.4 ng/dL
- Upper Limit of 6.0 ng/dL
- Assay time: 7 minutes
- CPT CODE: 84439

ORDERING INFORMATION	CONFIGURATION	CATALOG NUMBER
FastPack® IP Free Thyroxine (Free T4)	30 Tests	25000043
FastPack® IP Free T4 Calibration	Kit	25000016
FastPack® IP Free T4 Verification	Kit	25000031
FastPack® Control (US)	Kit	25000056
FastPack® Control (International)	Kit	25000065

THE AMERICAS

Sekisui Diagnostics, LLC
4 Hartwell Place
Lexington, MA 02421
Phone: 800 332 1042
Fax: 800 762 6311
Email: questions@sekisui-dx.com

INTERNATIONAL

Sekisui Diagnostics (UK) Limited
Liphook Way, Allington
Maidstone, Kent ME16 0LQ, UK
Phone: +44 1622 607800
Fax: +44 1622 607801
Email: info@sekisui-dx.com

SEKISUI
DIAGNOSTICS

Because every result matters™

www.sekisuidiagnostics.com