

The FastPack® IP PSA Immunoassay



FAST / SIMPLE / ACCURATE

The FastPack® IP PSA Immunoassay

The FastPack® IP PSA Immunoassay, for use on the FastPack® IP System, allows any testing facility to produce point-of-care results in 12 minutes with less than one minute of hands-on time.

- Enhanced quality of service for patients
- Ability to monitor and treat in one visit
- Improved patient compliance and outcome
- Laboratory accuracy and precision



SEKISUI
DIAGNOSTICS

Because every result matters™

The FastPack® IP PSA Immunoassay

Simple Process



1
*Dispense sample into the
FastPack® Injection Port*



2
*Place FastPack® into the
FastPack® IP Instrument*



3
*Press the **START** button*

The FastPack® IP PSA Assay Specifications:

- Quantitative measurement of Total PSA
- Sample: serum or plasma
- Lower limit: 0.04 ng/mL
- Upper limit: 50 ng/mL
- Assay time: 12 minutes
- CPT CODE: 84153

ORDERING INFORMATION	CONFIGURATION	CATALOG NUMBER
FastPack® IP Prostate Specific Antigen (PSA)	30 Tests	25000040
FastPack® IP PSA Calibration	Kit	25000002
FastPack® IP PSA Verification	Kit	25000017
FastPack® Control (US)	Kit	25000056
FastPack® PSA/Free PSA Control (International)	Kit	25000003

THE AMERICAS

Sekisui Diagnostics, LLC
4 Hartwell Place
Lexington, MA 02421
Phone: 800 332 1042
Fax: 800 762 6311
Email: questions@sekisui-dx.com

INTERNATIONAL

Sekisui Diagnostics (UK) Limited
Liphook Way, Allington
Maidstone, Kent ME16 0LQ, UK
Phone: +44 1622 607800
Fax: +44 1622 607801
Email: info@sekisui-dx.com

www.sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

Because every result matters™