

The FastPack® IP Vitamin D Immunoassay

FAST / SIMPLE / ACCURATE

The FastPack® IP Vitamin D Immunoassay

The FastPack® IP Vitamin D Immunoassay, for use on the FastPack® IP System, allows any testing facility to produce point-of-care results in 12 minutes with less than one minute of hands-on time.

- Enhanced quality of service for patients
- Ability to monitor and treat in one visit
- Improved patient compliance and outcome
- Laboratory accuracy and precision



SEKISUI
DIAGNOSTICS

Because every result matters™

The FastPack® IP Vitamin D Immunoassay

Simple Process

PRETREAT SAMPLE

- 1) *Pipette 100 µL of sample into prefilled vial*
- 2) *Recap and invert 3 times*



Dispense sample into the FastPack® Injection Port



Place FastPack® into the FastPack® IP Instrument



Press the START button

The FastPack® IP Vitamin D Assay Specifications:

- Quantitative measurement of 25-hydroxyvitamin D
- Lower Limit of 12.9 ng/dL
- Upper Limit of 150 ng/dL
- Assay time: 12 minutes
- CPT CODE: 82306

| ORDERING INFORMATION | CONFIGURATION | CATALOG NUMBER |
|---|---------------|----------------|
| FastPack® IP Vitamin D (Includes Calibrator, Controls, and Pre-Treatment Buffer) | 30 Tests | 25000068 |
| FastPack® IP Vitamin D Verification | Kit | 25000064 |

THE AMERICAS

Sekisui Diagnostics, LLC
 4 Hartwell Place
 Lexington, MA 02421
 Phone: 800 332 1042
 Fax: 800 762 6311
 Email: questions@sekisui-dx.com

INTERNATIONAL

Sekisui Diagnostics (UK) Limited
 Liphook Way, Allington
 Maidstone, Kent ME16 0LQ, UK
 Phone: +44 1622 607800
 Fax: +44 1622 607801
 Email: info@sekisui-dx.com

www.sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS
Because every result matters™