

RAPID DIAGNOSTICS

osom[®]

Strep A Test

FOR THE DETECTION OF GROUP A *STREPTOCOCCUS* ANTIGEN

www.osomtests.com



FEATURES:

- ✓ **ACCURATE**; 96% sensitivity; 98% specificity versus culture
- ✓ **FAST**; Results in 5 minutes
- ✓ **SAMPLE**; Throat swabs
- ✓ **QC INSIDE**; Includes 2 additional tests for external QC
- ✓ **TRAINING**; Online training and certification; www.osomtraining.com
- ✓ **CUSTOM CARE**; Technical support staffed with certified medical technologists



RAPID DIAGNOSTICS

www.osomtests.com

osom[®]

Strep A Test

Sekisui Diagnostics, LLC (AMERICAS)

4 Hartwell Place Lexington, MA 02421

Phone: 800 332 1042

Fax: 800 762 6311

Email: osomrapidtests@sekisui-dx.com

ISO 13485 Certified

Sekisui Diagnostics (UK) Limited (INTERNATIONAL)


Liphook Way, Allington, Maidstone, Kent ME16 0LQ, UK

Phone: +44 1622 607800 Fax: +44 1622 607801

Email: info@sekisui-dx.com

ISO 9001 Certified, ISO 13485 Certified

TEST PROCEDURE / RESULTS



DESCRIPTION	CONFIGURATION	CATALOG NUMBER	STORAGE	CPT CODE
OSOM [®] STREP A TEST	50 TESTS + 2 TESTS FOR QC	141 OR 141E (CE MARKED)	ROOM TEMP	87880QW
OSOM [®] STREP A TEST	20 TESTS + 2 TESTS FOR QC	141E-20 (CE MARKED)		

