



hCG Urine and Combo




Two simple tests you can depend on for the early detection of pregnancy


- Results in 5 minutes (serum), 3 minutes (urine), positive result may become visible sooner
- >99% sensitivity, specificity, accuracy
- 25 mIU/mL sensitivity in serum or urine
- Simple, 1-step procedure uses only 3 drops of sample
- Highly sensitive assay with built-in controls
- Highly controlled and monitored manufacturing environment brings lot-to-lot consistency
- Room temperature storage with 18-month shelf life from date of manufacture.



Combo Procedure:



Urine Procedure:



The performance of any rapid pregnancy test is dependent on sample collection and handling and the adherence to the Package Insert. Refer to our website at quidel.com for additional performance claims.

SL0803 | SS1008000EN00 (10/16)

QuickVue hCG Urine Test – 25 Test Kit: Cat. #20109
 QuickVue hCG Combo Test – 50 Test Kit: Cat. #20110
 hCG Control Set (Urine) – Cat. #00272
 hCG Control Set (Serum) – Cat. #00281

