

OPA/28 Disinfectant

McKesson OPA/28 disinfectant is a high-level, ready-to-use solution that provides a superior way to reprocess heat-sensitive, semi-critical medical devices that are unsuitable for sterilization. McKesson OPA/28 test strips verify the efficacy of McKesson OPA/28 solution. The user should test for the minimum effective concentration (MEC) prior to each use. McKesson OPA/28 test strips provides a fast disinfectant reaction time and user-friendly color comparison chart. McKesson OPA/28 disinfectant helps ensure a safe and healthy environment for patients and staff.

Features and Benefits

OPA/28 Disinfectant

- Excellent materials compatibility and biocompatibility
- Ready to use, no mixing or activating
- 10-minute manual disinfectant time at room temperature of 68°F (20°C)
- 5-minute AER disinfection time at 77°F (25°C)
- Not made with natural rubber latex

OPA/28 Disinfectant Test Strips

- 3-second dip time, read at 90 seconds
- Large, easy-to-read indicator pad
- Definitive color identification comparison chart on each bottle
- Not made with natural rubber latex



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
OPA/28, Test Strips	73-OPA28B	50 strips/bottle, 2 bottles/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.