

Multi-Enzymatic Cleanser

McKesson Multi-Enzymatic Cleanser contains a synergistic combination of multiple enzymes, which effectively cleans away blood, proteins, fat, carbohydrates, and synthetic liquids from surgical instruments, flexible and rigid scopes, and other general healthcare equipment. Its neutral pH, non-abrasive, non-foaming, 100% biodegradable formula is easy to use, and is specifically created to not harm metal, plastic, rubber, corrugated tubing on the equipment being cleaned. McKesson Multi-Enzymatic Cleanser is effective in both manual and automated reprocessing systems, as well as ultrasonic units and washer sterilizer/decontaminators.

Features and Benefits

- Rapidly dissolves fat, blood, mucus and other proteinaceous materials
- Liquid solution disperses instantly, covering entire surface of the instrument
- Non-corrosive formula with rinsing agents is safe for glass, plastic, rubber and metal surgical instruments and equipment
- Effective proteolytic detergent for the cleaning of flexible and rigid fiber optic scopes, delicate surgical instruments and other medical devices
- Low-foaming, neutral pH formula designed for pre-soaking and manual washing
- Phosphate-free, Chlorine-free, biodegradable



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Multi-Enzymatic Cleanser, 1 gal	53-28501	4/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.