

CAPIJECT® Safety Lancets

The Perfect Pairing of Ease and Precision.



- Ergonomically shaped for comfortable use
- Promotes confident operation and handling
- Single-hand, one-step operation
- Reduced training time
- No special staging, ready to use right out of the box
- Correct product selection is easier with color codes
- Top mounted activation button maximizes visibility for precise placement
- Specifically designed for single use only, preventing reuse
- Automatic retraction of lancet minimizes exposure of sharp and possible contamination

A legacy of proven design and trusted performance:

- Unmatched design
- Saves time
- Greater accuracy
- Promotes safety

Catalog Number	Lancet Color	Width/ Gauge	Depth	Blood Volume
200101	Pink	1.5mm	1.0mm	Low Blood flow
200102	Green	1.5mm	1.5mm	Medium Blood flow
200103	Blue	1.5mm	2.0mm	High Blood flow
200104	Orange	23G	2.25mm	Single Blood Drop
200105	Purple	28G	1.25mm	Single Blood Drop

Packaging:

- 200 Lancets per box
- 10 boxes per case-2,000 Lancets

These products contain no components made of natural rubber latex.

TerumoTMP.com/Capiject

TERUMO and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. CAPIJECT is a registered trademark of Terumo Medical Corporation, Somerset, NJ. TERUMO is registered trademark of Terumo Corporation, Tokyo, Japan ©2017 Terumo Medical Corporation. TMP-0263-01162017



TERUMO MEDICAL PRODUCTS