

PEEK-OPTIMA[®] HA Enhanced from Invibio[®] Biomaterial Solutions

Early, enhanced bone apposition

Product Summary

Invibio PEEK-OPTIMA[®] HA Enhanced provides the accepted high performance characteristics of PEEK-OPTIMA Natural, the leading spinal interbody device material, plus Hydroxyapatite (HA) to enhance bone apposition. This combination is an excellent solution for medical applications where early bone on-growth is required, such as spinal interbody fusion.

Product Information

- PEEK-OPTIMA[®] HA offers all the clinical advantages of PEEK-OPTIMA:
 - modulus similar to bone
 - reduced stress shielding
 - reduced subsidence
 - artifact free imaging
 - Plus the added benefit of enhanced bone apposition on all surfaces of the cage.
- PEEK-OPTIMA[®] HA Enhanced has been shown to improve bone apposition compared to unfilled PEEK-OPTIMA in a pre-clinical ovine model at 4 weeks and 12 weeks. These results do not necessarily correlate to human clinical outcomes.
- PEEK-OPTIMA[®] HA Enhanced makes hydroxyapatite available at all surfaces versus plasma spray coatings which typically have coverage on the endplates only.
- PEEK-OPTIMA[®] HA Enhanced has mechanical properties similar to PEEK-OPTIMA Natural and a modulus closer to bone than titanium, thus reducing the potential for stress shielding.
- PEEK-OPTIMA[®] HA Enhanced does not look different than PEEK-OPTIMA, but the hydroxyapatite particles, although naked to the eye, can be clearly seen on SEM images of the implant surface.

- PEEK-OPTIMA® HA Enhanced is radiolucent on X-rays and is artefact free on CT scans and MRI.
- PEEK-OPTIMA® HA Enhanced is an excellent choice for medical applications where early bone-on growth is required, such as spinal interbody fusion

Clinical Applications

The compounding of PEEK-OPTIMA HA Enhanced polymer with Hydroxyapatite (HA) enhances bone apposition, along with delivering the same high performance as PEEK-OPTIMA Natural even after repeated sterilization. It is an excellent choice for medical applications where early bone-on growth is required, such as spinal interbody fusion.



Europe and Asia Pacific
Invibio Ltd.
Technology Centre, Hillhouse International
Thornton Cleveleys, Lancashire FY5 4QD
United Kingdom
Tel: +44 (0) 1253 898000
FAX: +44 (0) 1253 898001

Americas
Invibio Inc.
300 Conshohocken State Road, Suite 120
West Conshohocken, PA 19428
USA
866-INVIBIO (468-4246)
Tel: +484 342 6004
FAX: +484 342 6005

For further information email us at info@invibio.com or please visit our website at:

► www.invibio.com

The information contained herein is believed to be an accurate description of the typical characteristics and/or uses of Invibio product(s). However, it is your ultimate responsibility to determine the performance, efficacy and safety of using Invibio product(s) for a specific application. Suggestions of uses should not be taken as inducements to infringe any particular patent or as a representation that the product is suitable for such uses. INVIBIO MAKES NO WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, A WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE OR OF INTELLECTUAL PROPERTY NON-INFRINGEMENT, INCLUDING, BUT NOT LIMITED TO PATENT NON-INFRINGEMENT, WHICH ARE EXPRESSLY DISCLAIMED, WHETHER EXPRESS OR IMPLIED, IN FACT OR BY LAW. FURTHER, INVIBIO MAKES NO WARRANTY TO YOUR CUSTOMERS OR AGENTS, AND HAS NOT AUTHORIZED ANYONE TO MAKE ANY REPRESENTATION OR WARRANTY OTHER THAN AS PROVIDED ABOVE.

INVIBIO SHALL IN NO EVENT BE LIABLE FOR ANY GENERAL, INDIRECT, SPECIAL, CONSEQUENTIAL, PUNITIVE, INCIDENTAL OR SIMILAR DAMAGES, INCLUDING WITHOUT LIMITATION, DAMAGES FOR HARM TO BUSINESS, LOST PROFITS OR LOST SAVINGS, EVEN IF INVIBIO HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, REGARDLESS TO THE FORM OF ACTION.

© 2014 Invibio Biomaterial Solutions.

Invibio® and the Invibio logo are registered trademarks of Invibio Ltd. PEEK-OPTIMA® is a registered trademark of Invibio Ltd.

S-MS-INV--E-0136-A