






CES[®]
CUTTING EDGE SPINE



STRAIGHT
LORDOTIC




TRUE
OBLIQUE
LORDOTIC



**TRUE OBLIQUE
LORDOTIC**

VOS HA
HA ENHANCED PEEK

Invio Biometrical Solutions' PEEK-OPTIMA® HA Enhanced is a material enhancement in spinal device technology. Hydroxyapatite (HA), a well-known osteo-conductive material, is fully integrated with strength, versatility, and performance advantages of its proven and popular predecessor. PEEK-OPTIMA® HA Enhanced offers a truly superior solution for bone apposition.



Rotate TOL	Heights (mm)	Widths (mm)	Lengths (mm)
Flat endplates, 5°TOL*	9, 10, 11	8	30
	10, 11, 12	9	30
	11, 12, 13	10	30
	12, 13, 14	11	30
	13, 14, 15, 16	12	30
Flat endplates, 8°TOL*	10, 11	8	30
	10, 11, 12	9	30
	11, 12, 13	10	30
	12, 13, 14	11	30
	13, 14, 15, 16	12	30

Straight TOL	Nose Styles	Heights (mm)	Widths (mm)	Lengths (mm)
Flat endplates, 5°TOL*	Standard, Wide, Narrow	8, 9, 10, 11, 12, 13, 14, 15, 16	8, 9, 10, 11	30

*True Oblique Lordotic