

SENTINEL_{plus}®

Enteral Feeding Pump

Become part of the evolution



Delivering the same
durability & reliability,
simply smaller.
Now that's a *plus*!

USA
MADE

ENTERAL FEEDING SYSTEMS FROM

ALCOR
SCIENTIFIC



DUAL DISPLAY

The bright LED display allows caregivers to simultaneously and conveniently view the infusion rate and volume delivered from a distance.



ADVANCED MEMORY

The pump operating history is retained until manually erased, even in the event of power failures or when pump is turned off.



DOSE CHECK

A record is kept of the accumulated amount of formula delivered over multiple feedings.



SIMPLE ALARMS

Audible and visual alarms are easily differentiated.



CLOVER-SHAPED ROTOR

The unique clover-shaped rotor allows for easy recognition of pump operation from a reasonable distance.



SENTINELplus® Advanced Technology

The SENTINELplus® Enteral Feeding Pump is a revised version of our popular SENTINEL® product. The lighter, more compact design brings the same simplicity and reliability of its predecessor. Contact us to learn more!

Order Information

Cat.	Description	Packaging
S-3000-IU	SENTINELplus®, Enteral Feeding Pump	1 each

Specifications

Mechanism	Rotary Peristaltic
Operating Range	5-295 mL/hour in 1 mL increments
Dose Limit	1-2000 mL in 1 mL increments 2000-9999 in 5 mL increments
Dose Accuracy	± 10% of selected flow rate with ALCOR sets
Occlusion Pressure	Does not exceed 15 psi
Battery	Battery backup 6 hours minimum at 125 mL/hr
Dimensions	15cm H x 20cm W x 9cm D
Weight	1.42kg
Power Options	All country-specific connectors available

www.alcorscientific.com

info@alcorscientific.com

401.737.3774

800.495.5270

113-09-500 2.3.17